




Phospho-DOK2 (Y299) Antibody

Product Code	CSB-PA000513
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60496
Immunogen	Synthesized peptide derived from Human Dok-2 around the phosphorylation site of Y299.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Monkey
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	DOK2; Docking protein 2; Downstream of tyrosine kinase 2; p56(dok-2)
Immunogen Species	Homo sapiens (Human)
Target Names	DOK2

Image



Western Blot analysis of COS7 cells using Phospho-Dok-2 (Y299) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.