





Phospho-FOXO1 (S319) Antibody

Product Code	CSB-PA000529
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12778
Immunogen	Synthesized peptide derived from Human FoxO1 around the phosphorylation site of S319.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	FOXO1; FKHR; FOXO1A; Forkhead box protein O1; Forkhead box protein O1A; Forkhead in rhabdomyosarcoma
Immunogen Species	Homo sapiens (Human)
Target Names	FOXO1

Image



Western Blot analysis of 293 cells using Phospho-FoxO1 (S319) Polyclonal Antibody



Western Blot analysis of NIH-3T3 cells using Phospho-FoxO1 (S319) Polyclonal Antibody