




Phospho-FOXO4 (S197) Antibody

Product Code	CSB-PA000530
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P98177
Immunogen	Synthesized peptide derived from Human FoxO4 around the phosphorylation site of S197.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	FOXO4; AFX; AFX1; MLLT7; Forkhead box protein O4; Fork head domain transcription factor AFX1
Immunogen Species	Homo sapiens (Human)
Target Names	FOXO4

Image



Western Blot analysis of COLO205 cells using Phospho-FoxO4 (S197) Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.