




Phospho-FOXO4 (S197) Antibody

Product Code	CSB-PA000530
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P98177
Immunogen	Synthesized peptide derived from Human FoxO4 around the phosphorylation site of S197.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	FOXO4; AFX; AFX1; MLLT7; Forkhead box protein O4; Fork head domain transcription factor AFX1
Immunogen Species	Homo sapiens (Human)
Target Names	FOXO4

Image



Western Blot analysis of COLO205 cells using Phospho-FoxO4 (S197) Polyclonal Antibody