




Phospho-GAB1 (Y627) Antibody

Product Code	CSB-PA000532
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13480
Immunogen	Synthesized peptide derived from Human Gab 1 around the phosphorylation site of Y627.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Monkey
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GAB1; GRB2-associated-binding protein 1; GRB2-associated binder 1; Growth factor receptor bound protein 2-associated protein 1
Immunogen Species	Homo sapiens (Human)
Target Names	GAB1

Image



Western Blot analysis of COS7 cells using Phospho-Gab 1 (Y627) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.