





# Phospho-HDAC4 (S632) Antibody

<b>Product Code</b>	CSB-PA000537
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P56524
<b>Immunogen</b>	Synthesized peptide derived from Human HDAC4 around the phosphorylation site of S632.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse, Rat
<b>Tested Applications</b>	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Isotype</b>	IgG
<b>Alias</b>	HDAC4; KIAA0288; Histone deacetylase 4; HD4
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	HDAC4

## Image



Western Blot analysis of 293 cells using Phospho-HDAC4 (S632) Polyclonal Antibody



Western blot analysis of HELA SH-SY5Y 3T3 lysis using Phospho-HDAC4 (S632) antibody.

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.