





Phospho-HDAC8 (S39) Antibody

Product Code	CSB-PA000539
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BY41
Immunogen	Synthesized peptide derived from Human HDAC8 around the phosphorylation site of S39.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	HDAC8; HDACL1; CDA07; Histone deacetylase 8; HD8
Immunogen Species	Homo sapiens (Human)
Target Names	HDAC8

Image



Western Blot analysis of VEC cells using Phospho-HDAC8 (S39) Polyclonal Antibody



Western Blot analysis of VEC 293 cells using Phospho-HDAC8 (S39) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.