





Phospho-Histone H2A.X (S139) Antibody

Product Code	CSB-PA000540
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16104
Immunogen	Synthesized peptide derived from Human Histone H2A.X around the phosphorylation site of S139.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	H2AFX; H2AX; Histone H2A.x; H2a/x
Immunogen Species	Homo sapiens (Human)
Target Names	H2AFX

Image



Western Blot analysis of 293 cells using Phospho-Histone H2A.X (S139) Polyclonal Antibody



Western Blot analysis of 4T1 VEC 293 cells using Phospho-Histone H2A.X (S139) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.