




Phospho-HSPB1 (S15) Antibody

Product Code	CSB-PA000544
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04792
Immunogen	Synthesized peptide derived from Human HSP27 around the phosphorylation site of S15.
Raised In	Rabbit
Species Reactivity	Human,Monkey
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	HSPB1; HSP27; HSP28; Heat shock protein beta-1; HspB1; 28 kDa heat shock protein; Estrogen-regulated 24 kDa protein; Heat shock 27 kDa protein; HSP 27; Stress-responsive protein 27; SRP27
Immunogen Species	Homo sapiens (Human)
Target Names	HSPB1

Image



Western Blot analysis of COS7 cells using Phospho-HSP27 (S15) Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.