




Phospho-ICAM1 (Y512) Antibody

Product Code	CSB-PA000547
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05362
Immunogen	Synthesized peptide derived from Human ICAM-1 around the phosphorylation site of Y512.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ICAM1; Intercellular adhesion molecule 1; ICAM-1; Major group rhinovirus receptor; CD antigen CD54
Immunogen Species	Homo sapiens (Human)
Target Names	ICAM1

Image



Western Blot analysis of A549(85) cells using Phospho-ICAM-1 (Y512) Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.