




Phospho-JAK2 (Y1007) Antibody

Product Code	CSB-PA000562
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60674
Immunogen	Synthesized peptide derived from Human JAK2 around the phosphorylation site of Y1007.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	JAK2; Tyrosine-protein kinase JAK2; Janus kinase 2; JAK-2
Immunogen Species	Homo sapiens (Human)
Target Names	JAK2

Image



Western Blot analysis of NIH-3T3 cells using Phospho-JAK2 (Y1007) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.