





Phospho-MAP2K2 (T394) Antibody

Product Code	CSB-PA000570
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P36507
Immunogen	Synthesized peptide derived from Human MEK-2 around the phosphorylation site of T394.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IP, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IP:2-5ug, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MAP2K2; MEK2; MKK2; PRKMK2; Dual specificity mitogen-activated protein kinase kinase 2; MAP kinase kinase 2; MAPKK 2; ERK activator kinase 2; MAPK/ERK kinase 2; MEK 2
Immunogen Species	Homo sapiens (Human)
Target Names	MAP2K2

Image



Western Blot analysis of VEC cells using Phospho-MEK-2 (T394) Polyclonal Antibody



Western Blot analysis of VEC 293 cells using Phospho-MEK-2 (T394) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.