




Phospho-MAP2K4 (S80) Antibody

Product Code	CSB-PA000572
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P45985
Immunogen	Synthesized peptide derived from Human MEK-4 around the phosphorylation site of S80.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MAP2K4; JNKK1; MEK4; MKK4; PRKMK4; SEK1; SERK1; SKK1; Dual specificity mitogen-activated protein kinase kinase 4; MAP kinase kinase 4; MAPKK 4; JNK-activating kinase 1; MAPK/ERK kinase 4; MEK 4; SAPK/ERK kinase 1; SEK1; Stress-activated pro
Immunogen Species	Homo sapiens (Human)
Target Names	MAP2K4

Image



Western Blot analysis of 293 cells using Phospho-MEK-4 (S80) Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.