




Phospho-ERBB2 (Y1221/Y1222) Antibody

Product Code	CSB-PA000577
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04626
Immunogen	Synthesized peptide derived from Human Neu around the phosphorylation site of Y1221/Y1222.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IP, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IP:2-5ug, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ERBB2; HER2; MLN19; NEU; NGL; Receptor tyrosine-protein kinase erbB-2; Metastatic lymph node gene 19 protein; MLN 19; Proto-oncogene Neu; Proto-oncogene c-ErbB-2; Tyrosine kinase-type cell surface receptor HER2; p185erbB2; CD antigen CD340
Immunogen Species	Homo sapiens (Human)
Target Names	ERBB2

Image



Western Blot analysis of HepG2 cells using Phospho-Neu (Y1221/Y1222) Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.