




Phospho-MAPK14 (T180) Antibody

Product Code	CSB-PA000593
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16539
Immunogen	Synthesized peptide derived from Human p38 around the phosphorylation site of T180.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat,Monkey
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MAPK14; CSBP; CSBP1; CSBP2; CSPB1; MXI2; SAPK2A; Mitogen-activated protein kinase 14; MAP kinase 14; MAPK 14; Cytokine suppressive anti-inflammatory drug-binding protein; CSAID-binding protein; CSBP; MAP kinase MXI2; MAX-interacting protein
Immunogen Species	Homo sapiens (Human)
Target Names	MAPK14

Image



Western Blot analysis of COS7 cells using Phospho-p38 (T180) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.