





Phospho-TP53 (S15) Antibody

Product Code	CSB-PA000595
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04637
Immunogen	Synthesized peptide derived from Human p53 around the phosphorylation site of S15.
Raised In	Rabbit
Species Reactivity	Human,Rat
Tested Applications	WB, IHC, IP, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IP:2-5ug, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TP53; P53; Cellular tumor antigen p53; Antigen NY-CO-13; Phosphoprotein p53; Tumor suppressor p53
Immunogen Species	Homo sapiens (Human)
Target Names	TP53

Image



Western Blot analysis of HT29 cells using Phospho-p53 (S15) Polyclonal Antibody



Western Blot analysis of 293T cells using Phospho-p53 (S15) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.