




# Phospho-TP53 (S20) Antibody

<b>Product Code</b>	CSB-PA000596
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P04637
<b>Immunogen</b>	Synthesized peptide derived from Human p53 around the phosphorylation site of S20.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat,Monkey
<b>Tested Applications</b>	WB,IHC,ELISA;WB:1/500-1/2000.IHC:1/100-1/300.ELISA:1/5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Isotype</b>	IgG
<b>Alias</b>	TP53; P53; Cellular tumor antigen p53; Antigen NY-CO-13; Phosphoprotein p53; Tumor suppressor p53
<b>Species</b>	Homo sapiens (Human)
<b>Target Names</b>	TP53

## Image



Western Blot analysis of COS7 cells using Phospho-p53 (S20) Polyclonal Antibody