





Phospho-TP73 (Y99) Antibody

Product Code	CSB-PA000606
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15350
Immunogen	Synthesized peptide derived from Human p73 around the phosphorylation site of Y99.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TP73; P73; Tumor protein p73; p53-like transcription factor; p53-related protein
Immunogen Species	Homo sapiens (Human)
Target Names	TP73

Image



Western Blot analysis of HeLa cells using Phospho-p73 (Y99) Polyclonal Antibody



Western Blot analysis of NIH-3T3 cells using Phospho-p73 (Y99) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.