




Phospho-PLCG1 (Y783) Antibody

Product Code	CSB-PA000612
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19174
Immunogen	Synthesized peptide derived from Human PLC γ 1 around the phosphorylation site of Y783.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat, Monkey
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PLCG1; PLC1; 1-phosphatidylinositol 4; 5-bisphosphate phosphodiesterase gamma-1; PLC-148; Phosphoinositide phospholipase C-gamma-1; Phospholipase C-II; PLC-II; Phospholipase C-gamma-1; PLC-gamma-1
Immunogen Species	Homo sapiens (Human)
Target Names	PLCG1

Image



Western Blot analysis of 293 cells using Phospho-PLC γ 1 (Y783) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.