





Phospho-STAT1 (Y701) Antibody

Product Code	CSB-PA000620
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42224
Immunogen	Synthesized peptide derived from Human Stat1 around the phosphorylation site of Y701.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat,Monkey
Tested Applications	WB, IHC, IP, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IP:2-5ug, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	STAT1; Signal transducer and activator of transcription 1-alpha/beta; Transcription factor ISGF-3 components p91/p84
Immunogen Species	Homo sapiens (Human)
Target Names	STAT1

Image



Western Blot analysis of COS7 cells using Phospho-Stat1 (Y701) Polyclonal Antibody



Western Blot analysis of MOUSE-BRAIN SH-SY5Y cells using Phospho-Stat1 (Y701) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.