





Phospho-STAT3 (S727) Antibody

Product Code	CSB-PA000621
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P40763
Immunogen	Synthesized peptide derived from Human Stat3 around the phosphorylation site of S727.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat,Monkey
Tested Applications	WB, IHC, IP, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IP:2-5ug, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	STAT3; APRF; Signal transducer and activator of transcription 3; Acute-phase response factor
Immunogen Species	Homo sapiens (Human)
Target Names	STAT3

Image



Western Blot analysis of COS7 cells using Phospho-Stat3 (S727) Polyclonal Antibody



Western Blot analysis of Jurkat cells using Phospho-Stat3 (S727) Polyclonal Antibody