





# Phospho-STAT3 (Y705) Antibody

<b>Product Code</b>	CSB-PA000622
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P40763
<b>Immunogen</b>	Synthesized peptide derived from Human Stat3 around the phosphorylation site of Y705.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse, Rat
<b>Tested Applications</b>	WB, IHC, IP, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IP:2-5ug, ELISA:1:20000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Isotype</b>	IgG
<b>Alias</b>	STAT3; APRF; Signal transducer and activator of transcription 3; Acute-phase response factor
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	STAT3

## Image



Western Blot analysis of KB cells using Phospho-Stat3 (Y705) Polyclonal Antibody



Western Blot analysis of KB MCF7 cells using Phospho-Stat3 (Y705) Polyclonal Antibody

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.