





Phospho-HDAC2 (S394) Antibody

Product Code	CSB-PA000633
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92769
Immunogen	Synthesized peptide derived from Human HDAC2 around the phosphorylation site of S394.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat,Monkey
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	HDAC2; Histone deacetylase 2; HD2
Immunogen Species	Homo sapiens (Human)
Target Names	HDAC2

Image



Western Blot analysis of 3T3 cells using Phospho-HDAC2 (S394) Polyclonal Antibody



Western Blot analysis of NIH-3T3 cells using Phospho-HDAC2 (S394) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.