





Phospho-BTK (Y551) Antibody

Product Code	CSB-PA000639
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q06187
Immunogen	Synthesized peptide derived from Human Btk around the phosphorylation site of Y551.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	BTK; AGMX1; ATK; BPK; Tyrosine-protein kinase BTK; Agammaglobulinaemia tyrosine kinase; ATK; B-cell progenitor kinase; BPK; Bruton tyrosine kinase
Immunogen Species	Homo sapiens (Human)
Target Names	BTK

Image



Western Blot analysis of HeLa cells using Phospho-Btk (Y551) Polyclonal Antibody



Western Blot analysis of HELA cells using Phospho-Btk (Y551) Polyclonal Antibody