




Phospho-RAD17 (S645) Antibody

Product Code	CSB-PA000641
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75943
Immunogen	Synthesized peptide derived from Human Rad17 around the phosphorylation site of S645.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	RAD17; R24L; Cell cycle checkpoint protein RAD17; hRad17; RF-C/activator 1 homolog
Immunogen Species	Homo sapiens (Human)
Target Names	RAD17

Image



Western Blot analysis of HeLa cells using Phospho-Rad17 (S645) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.