




# Phospho-RAD17 (S645) Antibody

<b>Product Code</b>	CSB-PA000641
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	O75943
<b>Immunogen</b>	Synthesized peptide derived from Human Rad17 around the phosphorylation site of S645.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse
<b>Tested Applications</b>	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Isotype</b>	IgG
<b>Alias</b>	RAD17; R24L; Cell cycle checkpoint protein RAD17; hRad17; RF-C/activator 1 homolog
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	RAD17

## Image



Western Blot analysis of HeLa cells using Phospho-Rad17 (S645) Polyclonal Antibody