





Phospho-OPRK1 (S369) Antibody

Product Code	CSB-PA000643
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41145
Immunogen	Synthesized peptide derived from Human KOR-1 around the phosphorylation site of S369.
Raised In	Rabbit
Species Reactivity	Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OPRK1; OPRK; Kappa-type opioid receptor; K-OR-1; KOR-1
Immunogen Species	Mus musculus (Mouse)
Target Names	OPRK1

Image



Western Blot analysis of Rat-kidney cells using Phospho-KOR-1 (S369) Polyclonal Antibody



Western Blot analysis of KB 22RV1 cells using Phospho-KOR-1 (S369) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.