





Phospho-PPARG (S112) Antibody

Product Code	CSB-PA000644
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P37231
Immunogen	Synthesized peptide derived from Human PPAR- γ around the phosphorylation site of S112.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PPARG; NR1C3; Peroxisome proliferator-activated receptor gamma; PPAR-gamma; Nuclear receptor subfamily 1 group C member 3
Immunogen Species	Homo sapiens (Human)
Target Names	PPARG

Image



Western Blot analysis of PC-3 cells using Phospho-PPAR- γ (S112) Polyclonal Antibody



Western Blot analysis of various cells using Phospho-PPAR- γ (S112) Polyclonal Antibody



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.