





Phospho-NOS3 (S1177) Antibody

Product Code	CSB-PA000659
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29474
Immunogen	Synthesized peptide derived from Human NOS3 around the phosphorylation site of S1177.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	NOS3; Nitric oxide synthase; endothelial; Constitutive NOS; cNOS; EC-NOS; Endothelial NOS; eNOS; NOS type III; NOSIII
Immunogen Species	Homo sapiens (Human)
Target Names	NOS3

Image



Western Blot analysis of HuvEc cells using Phospho-NOS3 (S1177) Polyclonal Antibody



Western Blot analysis of 22RV1 cells using Phospho-NOS3 (S1177) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.