





# Phospho-NOS3 (S1177) Antibody

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | CSB-PA000659  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.   |
| <b>Uniprot No.</b>         | P29474  |
| <b>Immunogen</b>           | Synthesized peptide derived from Human NOS3 around the phosphorylation site of S1177.                                 |
| <b>Raised In</b>           | Rabbit  |
| <b>Species Reactivity</b>  | Human,Mouse,Rat   |
| <b>Tested Applications</b> | WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:20000  |
| <b>Form</b>                | Liquid  |
| <b>Conjugate</b>           | Non-conjugated  |
| <b>Storage Buffer</b>      | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.   |
| <b>Purification Method</b> | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| <b>Isotype</b>             | IgG   |
| <b>Alias</b>               | NOS3; Nitric oxide synthase; endothelial; Constitutive NOS; cNOS; EC-NOS; Endothelial NOS; eNOS; NOS type III; NOSIII |
| <b>Immunogen Species</b>   | Homo sapiens (Human)  |
| <b>Target Names</b>        | NOS3  |

## Image



Western Blot analysis of HuvEc cells using Phospho-NOS3 (S1177) Polyclonal Antibody



Western Blot analysis of 22RV1 cells using  
Phospho-NOS3 (S1177) Polyclonal Antibody