




Phospho-RPS6KA5 (S376) Antibody

Product Code	CSB-PA000661
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75582
Immunogen	Synthesized peptide derived from Human MSK1 around the phosphorylation site of S376.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	RPS6KA5; MSK1; Ribosomal protein S6 kinase alpha-5; S6K-alpha-5; 90 kDa ribosomal protein S6 kinase 5; Nuclear mitogen- and stress-activated protein kinase 1; RSK-like protein kinase; RSKL
Immunogen Species	Homo sapiens (Human)
Target Names	RPS6KA5

Image



Western Blot analysis of HuvEc cells using Phospho-MSK1 (S376) Polyclonal Antibody