




Phospho-SYT1/SYT2 (T202/199) Antibody

Product Code	CSB-PA000678
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21579/Q8N9I0
Immunogen	Synthesized peptide derived from Human Synaptotagmin 1/2 around the phosphorylation site of T202/199.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	SYT1; SVP65; SYT; Synaptotagmin-1; Synaptotagmin I; SytI; p65; SYT2; Synaptotagmin-2; Synaptotagmin II; SytII
Immunogen Species	Homo sapiens (Human)
Target Names	SYT1/SYT2

Image



Western Blot analysis of 293 cells using Phospho-Synaptotagmin 1/2 (T202/199) Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.