




Phospho-LIMK1/LIMK2 (T508/505) Antibody

Product Code	CSB-PA000681
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P53667/P53671
Immunogen	Synthesized peptide derived from Human LIMK-1/2 around the phosphorylation site of T508/505.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	LIMK1; LIMK; LIM domain kinase 1; LIMK-1; LIMK2; LIM domain kinase 2; LIMK-2
Immunogen Species	Homo sapiens (Human)
Target Names	LIMK1/LIMK2

Image



Western Blot analysis of HeLa cells using Phospho-LIMK-1/2 (T508/505) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.