




Phospho-KSR1 (S392) Antibody

Product Code	CSB-PA000690
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IVT5
Immunogen	Synthesized peptide derived from Human Ksr-1 around the phosphorylation site of S392.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	KSR1; KSR; Kinase suppressor of Ras 1
Immunogen Species	Homo sapiens (Human)
Target Names	KSR1

Image



Western Blot analysis of HepG2 cells using Phospho-Ksr-1 (S392) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.