





Phospho-SYK (Y348) Antibody

Product Code	CSB-PA000697
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43405
Immunogen	Synthesized peptide derived from Human Syk around the phosphorylation site of Y348.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat,Monkey
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	SYK; Tyrosine-protein kinase SYK; Spleen tyrosine kinase; p72-Syk
Immunogen Species	Homo sapiens (Human)
Target Names	SYK

Image



Western Blot analysis of 3T3 cells using Phospho-Syk (Y348) Polyclonal Antibody



Western Blot analysis of NIH-3T3 cells using Phospho-Syk (Y348) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.