





Phospho-FOS (S362) Antibody

Product Code	CSB-PA000700
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01100
Immunogen	Synthesized peptide derived from Human c-Fos around the phosphorylation site of S362.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	FOS; G0S7; Proto-oncogene c-Fos; Cellular oncogene fos; G0/G1 switch regulatory protein 7
Immunogen Species	Homo sapiens (Human)
Target Names	FOS

Image



Western Blot analysis of 293T cells using Phospho-c-Fos (S362) Polyclonal Antibody



Western Blot analysis of 293T AD293 HELA HepG2 22RV1 cells using Phospho-c-Fos (S362) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.