





Phospho-PTK2 (Y407) Antibody

Product Code	CSB-PA000702
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q05397
Immunogen	Synthesized peptide derived from Human FAK around the phosphorylation site of Y407.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat,Monkey
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PTK2; FAK; FAK1; Focal adhesion kinase 1; FADK 1; Focal adhesion kinase-related nonkinase; FRNK; Protein phosphatase 1 regulatory subunit 71; PPP1R71; Protein-tyrosine kinase 2; p125FAK; pp125FAK
Immunogen Species	Homo sapiens (Human)
Target Names	PTK2

Image



Western Blot analysis of 293T cells using Phospho-FAK (Y407) Polyclonal Antibody



Western Blot analysis of 293T cells using
Phospho-FAK (Y407) Polyclonal Antibody