





Phospho-TP63 (S455) Antibody

Product Code	CSB-PA000705
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H3D4
Immunogen	Synthesized peptide derived from Human p63 around the phosphorylation site of S455.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TP63; KET; P63; P73H; P73L; TP73L; Tumor protein 63; p63; Chronic ulcerative stomatitis protein; CUSP; Keratinocyte transcription factor KET; Transformation-related protein 63; TP63; Tumor protein p73-like; p73L; p40; p51
Immunogen Species	Homo sapiens (Human)
Target Names	TP63

Image



Western Blot analysis of 293T cells using Phospho-p63 (S455) Polyclonal Antibody



Western Blot analysis of 293T cells using Phospho-p63 (S455) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.