





Phospho-PIK3R1 (Y607) Antibody

Product Code	CSB-PA000712
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P27986
Immunogen	Synthesized peptide derived from Human PI 3-kinase p85 α around the phosphorylation site of Y607.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PIK3R1; GRB1; Phosphatidylinositol 3-kinase regulatory subunit alpha; PI3-kinase regulatory subunit alpha; PI3K regulatory subunit alpha; PtdIns-3-kinase regulatory subunit alpha; Phosphatidylinositol 3-kinase 85 kDa regulatory subunit alph
Immunogen Species	Homo sapiens (Human)
Target Names	PIK3R1

Image



Western Blot analysis of KB cells using Phospho-PI 3-kinase p85 α (Y607) Polyclonal Antibody



Western Blot analysis of KB cells using Phospho-PI 3-kinase p85 α (Y607) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.