




Phospho-DDR1 (Y513) Antibody

Product Code	CSB-PA000715
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q08345
Immunogen	Synthesized peptide derived from Human DDR1 around the phosphorylation site of Y513.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	DDR1; CAK; EDDR1; NEP; NTRK4; PTK3A; RTK6; TRKE; Epithelial discoidin domain-containing receptor 1; Epithelial discoidin domain receptor 1; CD167 antigen-like family member A; Cell adhesion kinase; Discoidin receptor tyrosine kinase; HGK2;
Immunogen Species	Homo sapiens (Human)
Target Names	DDR1

Image



Western Blot analysis of HepG2 cells using Phospho-DDR1 (Y513) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.