




Phospho-APP (T743) Antibody

Product Code	CSB-PA000724
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05067
Immunogen	Synthesized peptide derived from Human Amyloid-β around the phosphorylation site of T743.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	APP; A4; AD1; Amyloid beta A4 protein; ABPP; APPI; APP; Alzheimer disease amyloid protein; Cerebral vascular amyloid peptide; CVAP; PreA4; Protease nexin-II; PN-II
Immunogen Species	Homo sapiens (Human)
Target Names	APP

Image



Western Blot analysis of HeLa cells using Phospho-Amyloid-β (T743) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.