





# Phospho-HSP90AB1 (S254) Antibody

<b>Product Code</b>	CSB-PA000726
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P08238
<b>Immunogen</b>	Synthesized peptide derived from Human HSP90β around the phosphorylation site of S254.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse, Rat, Monkey
<b>Tested Applications</b>	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:40000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Isotype</b>	IgG
<b>Alias</b>	HSP90AB1; HSP90B; HSPC2; HSPCB; Heat shock protein HSP 90-beta; HSP 90; Heat shock 84 kDa; HSP 84; HSP84
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	HSP90AB1

## Image



Western Blot analysis of HepG2 cells using Phospho-HSP90β (S254) Polyclonal Antibody



Western Blot analysis of K562 cells using  
Phospho-HSP90β (S254) Polyclonal Antibody