





Phospho-HSP90AB1 (S254) Antibody

Product Code	CSB-PA000726
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08238
Immunogen	Synthesized peptide derived from Human HSP90β around the phosphorylation site of S254.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat, Monkey
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	HSP90AB1; HSP90B; HSPC2; HSPCB; Heat shock protein HSP 90-beta; HSP 90; Heat shock 84 kDa; HSP 84; HSP84
Immunogen Species	Homo sapiens (Human)
Target Names	HSP90AB1

Image



Western Blot analysis of HepG2 cells using Phospho-HSP90β (S254) Polyclonal Antibody



Western Blot analysis of K562 cells using Phospho-HSP90 β (S254) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.