




Phospho-NFKB1 (S927) Antibody

Product Code	CSB-PA000730
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19838
Immunogen	Synthesized peptide derived from Human NFκB-p105 around the phosphorylation site of S927.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	WB,IHC,IF,ELISA;WB:1/500-1/2000.IHC:1/100-1/300.IF:1/200-1/1000.ELISA:1/20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	NFKB1; Nuclear factor NF-kappa-B p105 subunit; DNA-binding factor KBF1; EBP-1; Nuclear factor of kappa light polypeptide gene enhancer in B-cells 1
Species	Homo sapiens (Human)
Target Names	NFKB1

Image



Western Blot analysis of HeLa cells using Phospho-NFκB-p105 (S927) Polyclonal Antibody