





Phospho-TH (S19) Antibody

Product Code	CSB-PA000731
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07101
Immunogen	Synthesized peptide derived from Human TH around the phosphorylation site of S19.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TH; TYH; Tyrosine 3-monooxygenase; Tyrosine 3-hydroxylase; TH
Immunogen Species	Homo sapiens (Human)
Target Names	TH

Image



Western Blot analysis of 293 cells using Phospho-TH (S19) Polyclonal Antibody



Western Blot analysis of VEC 293 A549 cells using Phospho-TH (S19) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.