




Phospho-ESPL1 (S801) Antibody

Product Code	CSB-PA000737
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14674
Immunogen	Synthesized peptide derived from Human Separase around the phosphorylation site of S801.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ESPL1; ESP1; KIAA0165; Separin; Caspase-like protein ESPL1; Extra spindle poles-like 1 protein; Separase
Immunogen Species	Homo sapiens (Human)
Target Names	ESPL1

Image



Western Blot analysis of A431 cells using Phospho-Separase (S801) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.