





Phospho-CTNND1 (Y228) Antibody

Product Code	CSB-PA000741
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60716
Immunogen	Synthesized peptide derived from Human p120 around the phosphorylation site of Y228.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CTNND1; KIAA0384; Catenin delta-1; Cadherin-associated Src substrate; CAS; p120 catenin; p120(ctn); p120(cas)
Immunogen Species	Homo sapiens (Human)
Target Names	CTNND1

Image



Western Blot analysis of 293T cells using Phospho-p120 (Y228) Polyclonal Antibody



Western Blot analysis of 293T cells using
Phospho-p120 (Y228) Polyclonal Antibody