




Phospho-HNF4A (S313) Antibody

Product Code	CSB-PA000743
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41235
Immunogen	Synthesized peptide derived from Human HNF4- α around the phosphorylation site of S313.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	HNF4A; HNF4; NR2A1; TCF14; Hepatocyte nuclear factor 4-alpha; HNF-4-alpha; Nuclear receptor subfamily 2 group A member 1; Transcription factor 14; TCF-14; Transcription factor HNF-4
Immunogen Species	Homo sapiens (Human)
Target Names	HNF4A

Image



Western Blot analysis of HuvEc cells using Phospho-HNF4- α (S313) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.