





Phospho-OPRM1 (S375) Antibody

Product Code	CSB-PA000744
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35372
Immunogen	Synthesized peptide derived from Human MOR-1 around the phosphorylation site of S375.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OPRM1; MOR1; Mu-type opioid receptor; M-OR-1; MOR-1; Mu opiate receptor; Mu opioid receptor; MOP; hMOP
Immunogen Species	Homo sapiens (Human)
Target Names	OPRM1

Image



Western Blot analysis of Lovo cells using Phospho-MOR-1 (S375) Polyclonal Antibody



Western Blot analysis of LOVO cells using Phospho-MOR-1 (S375) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.