





# Phospho-FLT1 (Y1048) Antibody

<b>Product Code</b>	CSB-PA000747
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P17948
<b>Immunogen</b>	Synthesized peptide derived from Human Flt-1 around the phosphorylation site of Y1048.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	WB, IHC, ELISA; IHC:1:100-1:300, ELISA:1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Isotype</b>	IgG
<b>Alias</b>	FLT1; FLT; FRT; VEGFR1; Vascular endothelial growth factor receptor 1; VEGFR-1; Fms-like tyrosine kinase 1; FLT-1; Tyrosine-protein kinase FRT; Tyrosine-protein kinase receptor FLT; FLT; Vascular permeability factor receptor
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	FLT1

## Image



Western Blot analysis of COLO 823 293T-UV cells using Phospho-Flt-1 (Y1048) Polyclonal Antibody



Western blot analysis of 293T lysis using Phospho-Flt-1 (Y1048) antibody.

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.